

1 Title: The Detection of SARS-CoV-2 using the Cepheid Xpert Xpress SARS-CoV-2 and Roche  
2 cobas SARS-CoV-2 Assays

3

4 Authors: Angelica Moran<sup>a\*</sup>, Kathleen G. Beavis<sup>a\*</sup>, Scott M. Matushek<sup>b</sup>, Carol Ciaglia<sup>b</sup>, Nina  
5 Francois<sup>b</sup>, Vera Tesic<sup>a</sup>, Nedra Love<sup>b</sup>

6 Author Affiliation:

7 <sup>a</sup>Department of Pathology, The University of Chicago, 5841 South Maryland Avenue, Chicago,  
8 Illinois 60637

9 <sup>b</sup>Clinical Microbiology Laboratory, University of Chicago Medicine, 5841 South Maryland  
10 Avenue, Chicago, Illinois 60637

11 Angelica Moran and Kathleen G. Beavis *contributed equally to this work. Author order was*  
12 *determined in order of increasing seniority.*

13

14 Corresponding Author: Kathleen G. Beavis, kbeavis@uchicago.edu

15

16

17

18

19

20

21

22

23 SARS-CoV-2, a novel coronavirus responsible for a December 2019 outbreak in Wuhan, China,  
24 causes a syndrome characterized by fever, cough, and dyspnea progressing to acute respiratory  
25 distress syndrome (1). SARS-CoV-2 quickly spread to other countries, with the new coronavirus  
26 2019 (COVID-19) disease declared a pandemic in March 2020 (2-4). Rapid testing for SARS-  
27 CoV-2 is important for epidemiological tracking and institution of quarantine procedures (5).  
28 The clinical description of COVID-19 continues to evolve; with transmission by asymptomatic  
29 individuals reported (6-8), widespread testing is necessary.

30

31 Multiple RT-PCR assays have received Emergency Use Authorization from the US Food and  
32 Drug Administration. The Roche cobas SARS-CoV-2 assay is a qualitative test that detects the  
33 SARS-CoV-2-specific ORF1 and part of the E-gene conserved across sarbecoviruses, including  
34 SARS-CoV-2 (9). The Cepheid Xpert Xpress SARS-CoV-2 assay also detects the pan-  
35 sarbecovirus E-gene but detects the N2 region of the N-gene as its SARS-CoV-2-specific target  
36 (10). This report compares results from specimens tested with both assays.

37

38 Eight nasal and 95 nasopharyngeal specimens were collected from inpatients and ambulatory  
39 patients at the University of Chicago. Samples were tested using the Roche cobas SARS-CoV-2  
40 assay on the cobas 6800 system (Roche Molecular Systems, Branchburg, NJ) and the Cepheid  
41 Xpert Xpress SARS-CoV-2 assay on the GeneXpert system (Cepheid, Sunnyvale, CA). Of these  
42 103 specimens, 42 tested positive and 60 tested negative with both systems for agreement of  
43 99%. Testing was repeated on the single specimen with discrepant results. For this specimen, the  
44 Roche assay was repeatedly negative for SARS-CoV-2. The initial Cepheid assay result was  
45 positive for SARS-CoV-2, though the cycle threshold (Ct) value for detection of the E-gene was

46 0.0 (negative) and 42.0 (low positive) for the N-gene. Repeat Cepheid testing was negative for  
47 both targets. These results suggest that SARS-CoV-2 was present at a very low concentration,  
48 near the detection limit of the Cepheid assay.

49

50 For the 42 positive samples, Ct values for the E-gene ranged from 15.05 to 39.75 (mean 26.35,  
51 SD 6.69) for the Roche assay and 13.6 to 38.2 (mean 24.8, SD 7.1) for the Cepheid assay. Using  
52 Bland-Altman analysis to assess agreement, Ct values were lower on the Cepheid assay in 37 of  
53 42 samples (mean difference -1.57, 95% limits of agreement -5.34, 2.20). This could be due to  
54 differences in primer sequences against the E-gene, reagents, or amplification conditions.

55

56 Limitations of this study include small sample size of SARS-CoV-2-positive specimens, as  
57 testing was limited to patients within our institution. The assays also detect different SARS-  
58 CoV-2-specific genes, which could lead to false-negative results if a mutation prevents primer  
59 binding. The Cepheid assay is a 45-minute random-access assay, with throughput dependent on  
60 number of instrument slots. The Roche platform is batch-based, accommodating 90 samples/run  
61 every 90 minutes. As each run requires up to three hours and 45 minutes, throughput is  
62 approximately 1 result/minute. Overall, the Cepheid Xpert Xpress SARS-CoV-2 and Roche  
63 cobas SARS-CoV-2 assays show excellent agreement (>99%), and their combined usage can be  
64 tailored to maximize SARS-CoV-2 testing.

65

66

67

68

## 69 References.

- 70 1. Chen N, Zhou M, Dong X, Qu J, Gong F, Han Y, Qiu Y, Wang J, Liu Y, Wei Y, Xia J,  
71 Yu T, Zhang X, Zhang L. 2020. Epidemiological and clinical characteristics of 99 cases  
72 of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. *Lancet*  
73 395:507-513.
- 74 2. Ki M, Task Force for -nCo V. 2020. Epidemiologic characteristics of early cases with  
75 2019 novel coronavirus (2019-nCoV) disease in Korea. *Epidemiol Health* 42:e2020007.
- 76 3. Spiteri G, Fielding J, Diercke M, Campese C, Enouf V, Gaynard A, Bella A,  
77 Sognamiglio P, Sierra Moros MJ, Riutort AN, Demina YV, Mahieu R, Broas M, Bengner  
78 M, Buda S, Schilling J, Filleul L, Lepoutre A, Saura C, Mailles A, Levy-Bruhl D,  
79 Coignard B, Bernard-Stoecklin S, Behillil S, van der Werf S, Valette M, Lina B,  
80 Riccardo F, Nicastrì E, Casas I, Larrauri A, Salom Castell M, Pozo F, Maksyutov RA,  
81 Martin C, Van Ranst M, Bossuyt N, Siira L, Sane J, Tegmark-Wisell K, Palmerus M,  
82 Broberg EK, Beaute J, Jorgensen P, Bundle N, Pereyaslov D, Adlhoch C, Pukkila J,  
83 Pebody R, Olsen S, et al. 2020. First cases of coronavirus disease 2019 (COVID-19) in  
84 the WHO European Region, 24 January to 21 February 2020. *Euro Surveill* 25. doi:  
85 10.2807/1560-7917.ES.2020.25.9.2000178.
- 86 4. Holshue ML, DeBolt C, Lindquist S, Lofy KH, Wiesman J, Bruce H, Spitters C, Ericson  
87 K, Wilkerson S, Tural A, Diaz G, Cohn A, Fox L, Patel A, Gerber SI, Kim L, Tong S, Lu  
88 X, Lindstrom S, Pallansch MA, Weldon WC, Biggs HM, Uyeki TM, Pillai SK,  
89 Washington State -nCoV VCIT. 2020. First Case of 2019 Novel Coronavirus in the United  
90 States. *N Engl J Med* 382:929-936.
- 91 5. Deng SQ, Peng HJ. 2020. Characteristics of and Public Health Responses to the  
92 Coronavirus Disease 2019 Outbreak in China. *J Clin Med* 9. doi: 10.3390/jcm9020575
- 93 6. Bai Y, Yao L, Wei T, Tian F, Jin DY, Chen L, Wang M. 2020. Presumed Asymptomatic  
94 Carrier Transmission of COVID-19. *JAMA* doi:10.1001/jama.2020.2565.
- 95 7. Rothe C, Schunk M, Sothmann P, Bretzel G, Froeschl G, Wallrauch C, Zimmer T, Thiel  
96 V, Janke C, Guggemos W, Seilmaier M, Drosten C, Vollmar P, Zwirgmaier K, Zange S,  
97 Wolfel R, Hoelscher M. 2020. Transmission of 2019-nCoV Infection from an  
98 Asymptomatic Contact in Germany. *N Engl J Med* 382:970-971.
- 99 8. Borges do Nascimento IJ, Cacic N, Abdulazeem HM, von Groote TC, Jayarajah U,  
100 Weerasekara I, Esfahani MA, Civile VT, Marusic A, Jeroncic A, Carvas Junior N, Pericic  
101 TP, Zakarija-Grkovic I, Meirelles Guimaraes SM, Luigi Bragazzi N, Bjorklund M, Sofi-  
102 Mahmudi A, Altujjar M, Tian M, Arcani DMC, O'Mathuna DP, Marcolino MS. 2020.  
103 Novel Coronavirus Infection (COVID-19) in Humans: A Scoping Review and Meta-  
104 Analysis. *J Clin Med* 9. doi: 10.3390/jcm9040941
- 105 9. Roche Molecular Systems. 2020. cobas SARS-CoV-2: package insert. U.S. Food and  
106 Drug Administration website. <https://www.fda.gov/media/136049/download> Accessed 8  
107 April 2020.
- 108 10. Cepheid. 2020. Xpert Xpress SARS-CoV-2 Test: package insert. U.S. Food and Drug  
109 Administration website. <https://www.fda.gov/media/136314/download> Accessed 8 April  
110 2020.

111

112